



TECHNISCHE  
UNIVERSITÄT  
WIEN

DISSERTATION

# Cool Science

ausgeführt am Atominstutut



der Technische Universität Wien  
Fakultät für Physik

unter der Anleitung von  
**Univ.Prof. Dipl.-Ing. Dr.techn. Gorge Hammond**  
und  
**Projektass. Dr.rer.nat Rodney MacKay MSc.**  
**Projektass. Dr.techn. Dr.techn. Dr.techn. Dipl.-Ing.**  
**Samantha Carter**

durch

**Daniel Jackson**

Matrikelnummer: 9-18-27-15-21-36  
Stadionallee 2  
1020 Wien

Wien, am 08.04.2020

*“The Setesh guard’s nose drips.”*  
TEAL’C

# Contents

<b>1 Cicero Word Generator</b>	<b>1</b>
1.1 Installation of National Instruments drivers . . . . .	1
1.2 Installation of National Instruments Cards . . . . .	1
1.3 Configuring Atticus . . . . .	1
1.3.1 Configure hardware timing / synchronization . . . . .	2
1.4 Configuration and Basic Usage of Cicero . . . . .	3
1.5 Saving of Settings and Sequences . . . . .	3
1.6 Sequence length limit . . . . .	3
<b>2 Electron beam setup</b>	<b>4</b>
2.1 Charatarization of a working CRT . . . . .	4
2.2 High Voltage Power Supply HVPS . . . . .	7
2.2.1 Ripple measurement . . . . .	7
2.3 CRT wiring . . . . .	8
2.4 Heater . . . . .	9
<b>Todo list</b>	<b>11</b>
<b>References</b>	<b>12</b>

# <sup>1</sup> 1 Cicero Word Generator

<sup>2</sup> This chapter describes the installation and initial setup of Cicero Word Generator<sup>[1]</sup>  
<sup>3</sup> on a PC running Windows 10 with analog and digital cards from National Instruments  
<sup>4</sup> (NI). The code is freely available on Github<sup>[2]</sup>. This chapter contains only differences,  
<sup>5</sup> problems, and possible solutions encountered when Cicero was installed for the PC  
<sup>6</sup> ‘Fritz Fantom’ which will be used for the QuaK experiment. It is therefore advised  
<sup>7</sup> to use the technical and user manual<sup>[3]</sup> in conjunction. The titles in this chapter and  
<sup>8</sup> font style with Courier and Boldface was mirrored to fit the manual.

## <sup>9</sup> 1.1 Installation of National Instruments drivers

<sup>10</sup> Before setting up the Cicero Word Generator, it is necessary to install the newest  
<sup>11</sup> .NET Framework<sup>[4]</sup> from Microsoft. For the first installation of NI drivers, NI-DAQmx  
<sup>12</sup> (version 9.3), NI-VISA (newest version), and NI-4888.2 (newest version) should be  
<sup>13</sup> downloaded from the National Instruments website<sup>[5]</sup>. When installing the NI drivers  
<sup>14</sup> it is possible to get an ‘Runtime Error!’. In this case it is necessary to set the Regional  
<sup>15</sup> format settings of Windows 10 to ‘English (United States)’<sup>[6]</sup>.

## <sup>16</sup> 1.2 Installation of National Instruments Cards

<sup>17</sup> After installation of the necessary drivers, the physical cards can be inserted into the  
<sup>18</sup> PCIe slots on the motherboard. On ‘Fritz Fantom’ the digital card (NI PCIe-6537B)  
<sup>19</sup> was installed in PCIe bus 3 while the analog cards (NI PCIe-6738) were installed in  
<sup>20</sup> PCIe bus 4 and 5.

## <sup>21</sup> 1.3 Configuring Atticus

<sup>22</sup> After installation of the NI cards, Atticus should be launched for the first time and  
<sup>23</sup> closed without changing any settings. After this, the NI-DAQmx drivers should be

updated to the newest version. If version 9.3 was not used when launching Atticus in this step, it could result in an error. After this, “Configuring Atticus” on the user manual can be followed. The **Server Name** was set to ‘Fritz\_Phantom’. **Dev1** to **Dev3** were set in the same ascending order as the physical installation on the motherboard.

1

2

change server n  
in lab? Fantom  
4 Phantom

5

### 1.3.1 Configure hardware timing / synchronization

For synchronization, a **Shared Sample Clock** was used with **Dev1** being the master card. The settings are summarized in table 1.1 and table 1.2. For **Dev3** ‘SampleClockExternalSource’ should be set to ‘/Dev3/RTSI7’. The ‘SampleClockRate’ is set to 350 kHz since this is the fastest rate with all 32 analog channels active. It is possible to raise this to 1 MHz by only using 8 channels (1 channel per bank).

6

7

8

9

10

11

**Table 1.1:** Settings for **Dev1**.

Setting	Value
MasterTimebaseSource	
MySampleClockSource	DerivedFromMaster
SampleClockRate	350000
UsingVariablenamebase	False
SoftTriggerLast	True
StartTriggerType	SoftwareTrigger

**Table 1.2:** Settings for **Dev2**.

Setting	Value
MasterTimebaseSource	
MySampleClockSource	External
SampleClockExternalSource	/Dev2/RTSI7
SampleClockRate	350000
UsingVariablenamebase	False
SoftTriggerLast	False
StartTriggerType	SoftwareTrigger

## **1.4 Configuration and Basic Usage of Cicero**

<sup>2</sup> After setting up the Atticus server, Cicero can be configured. In step 3.c. it is necessary  
<sup>3</sup> to write the full IP address and not ‘localhost’. Once step 6 is finished, Cicero should  
<sup>4</sup> run without any problems.

## **1.5 Saving of Settings and Sequences**

<sup>6</sup> The ‘SettingsData’ of the Server Atticus are saved in C:\Users\confetti\Documents  
<sup>7</sup> \Cicero\_Atticus\Cicero\SettingsData while the ‘SequenceData’ of Cicero are saved in  
<sup>8</sup> C:\Users\confetti\Documents\Cicero\_Atticus\Cicero\SequenceData.

## **1.6 Sequence length limit**

<sup>10</sup> The duration of a sequence is limited to  $2^{32}/(16 * 32 * 350 \text{ kHz}) = 23.967 \text{ s}$  coming  
<sup>11</sup> from a 32-bit application, 16 bit per channel, 32 channels in a NI PCIe-6738 card, and  
<sup>12</sup> 350 kHz clock rate.

## 2 Electron beam setup

1

### 2.1 Charatarization of a working CRT

2

HAMEG HM507 oscilloscopes [7] were used for testing purposes. These contain a D14-363GY/123[8] CRT hereinafter abbreviated as ‘D14’, ‘tube’, or ‘CRT’. Although the HM507 has only a bandwidth of 0 MHz to 50 MHz, which is not sufficient for the hyperfine splitting frequency of 461.7 MHz of  $^{39}\text{K}$ , it was used nevertheless because of its simple construction and availability. A schematic view of the device is shown in fig. 2.1 with the back pin arrangement in fig. 2.2.

3

4

5

<http://www.to>

6

7

8

The voltages and currents of the necessary pins to drive the CRT were measured using a 2.5 kV probe with an attenuation ratio of and are summarized in table 2.1. It was not possible to measure pin g3 directly. Therefore a HVPS (section 2.2) was used to set a voltage and the beam diameter was observed. The best focus was achieved with the voltage mentioned in the table. The voltage offset of x-, and y-plates was not possible to measure directly, since it varies with time to draw the necessary image on the phosphor screen. The given values in table 2.1 are the mean of the minimum and maximum measured voltage. The deflection coefficient is summarized in table 2.2.

9

10model number

11:100 or 100:1

12

13current?


14

15

16


## 2 Electron beam setup

---



**Figure 2.1:** Electrode configuration (from [8])

how to cite figure



**Figure 2.2:** Pin arrangement, bottom view (from [8])

how to cite figure

## 2 Electron beam setup

---

**Table 2.1:** D14-363GY/123 CRT pin measurements

current empty or '-' symbol

number	pin	voltage/V	current/ $\mu$ A
1	f	$-1.99 \times 10^3$	$86.6 \times 10^3$
2	k	-2.00	-7.6
3	g1	-2.03	0
4	g3	$-1.813 \times 10^3$	
5	i.c.	71.7	0.1
6	g5	64.0	7.2
7	y2	78	
9	y1	78	
11	x2	96	-
12	g2, g4	71.0	0
13	x1	96	-
14	f	$-1.97 \times 10^3$	$-86.2 \times 10^3$

**Table 2.2:** D14-363GY/123 deflection coefficient (from [8])

how to cite source


horizontal	$M_x$	19 V/cm
vertical	$M_y$	11.5 V/cm

## <sup>1</sup> 2.2 High Voltage Power Supply HVPS

<sup>2</sup> To produce high dc voltages to drive the CRT, four HCP 14-6500 power supplies<sup>[9]</sup>  
<sup>3</sup> were used. They were named ‘HVPS 1’ to ‘HVPS 4’ and can provide up to  $\pm 6.5$  kV and  
<sup>4</sup> 2 mA. To connect the output to the CRT pins, BNC cables were refitted with a save  
<sup>5</sup> high voltage (SHV) connector on one side while on the other end the BNC connector  
<sup>6</sup> was kept (fig. 2.3). A 6 kV probe was used to obtain the breakdown voltage, which is  
<sup>7</sup> around 3 kV caused by the coaxial cable which was not built do sustain high voltages.

find name of bi  
yellow probe

somewhere 2.5-  
find exact value



**Figure 2.3:** Coaxial cable with SHV and BNC connector.


### <sup>8</sup> 2.2.1 Ripple measurement

<sup>9</sup> Each power supply was measured for its ripple with a set voltage of 2 kV. A 2.5 kV  
<sup>10</sup> probe (attenuation ratio )was connected to an oscilloscope set to ac coupling with a  
<sup>11</sup> timescale of 1 ms. To get the electronic noise of the oscilloscope itself, the probe was  
<sup>12</sup> shorted and the noise measured. A picture of a measurement is shown in fig. 2.4 with  
<sup>13</sup> the values summarized in table 2.3.

100:1 or 1:100

**Table 2.3:** HVPS ripple

device	ripple/mV
short	116
HVPS 1	136
HVPS 2	138
HVPS 3	194
HVPS 4	204




**Figure 2.4:** Measurement of HVPS ripple.

## 2.3 CRT wiring

A schematic of the supplied power is shown in fig. 2.5. A small ac or dc voltage is necessary to drive the heater filament f. This part of the setup is explained in section 2.4.


1  
2  
3  
4



**Figure 2.5:** Schematics of supplying CRT pins with power.

## <sup>1</sup> 2.4 Heater

- <sup>2</sup> The heater provides an adjustable ac voltage, which is used to regulate the temperature  
<sup>3</sup> of the cathode. In the cold state, the heater filament has a an electrical resistitance  
<sup>4</sup> of approximately  $15\Omega$ , when the filament is hot, this value rises to  $90\Omega$ . The normal  
<sup>5</sup> heater voltage for the D14-363GY/123 during operation is 6.0 V to 6.6 V according to  
<sup>6</sup> [8]. Our ac-power supply (figure 2.6 shows its circuit diagram) consists of an isolation



**Figure 2.6:** Circuit diagram of filament power supply.

check if really  $14\Omega$  or if it event exists

transformer (from grid voltage to 12 V), its primary and secondary circuits are isolated up to 4 kV [10]. The power supply has two banana plug sockets to connect to the heater filament. It is connected to the transformer in series with a  $100\Omega$  potentiometer. Using the full resistance, there is a voltage of approximately 5.7 V applied to the heater filament, by lowering the resistance this value can goes up to nearly the full voltage of the transformer. The current running through the filament is measured with an integrated amperemeter [11] that measures currents up to two 2 A with mA accuracy.

At the beginning of operation it is recommendable to set the maximum resistance and slowly increase the current to the desired value once the filament is heated up. As the resistance of the cold filament is significantly lower, high onset currents could otherwise damage it.

# <sup>1</sup> Todo list

<sup>2</sup>	■ namechange? . . . . .	2
<sup>3</sup>	■ <a href="http://www.tobiastiecke.nl/archive/PotassiumProperties.pdf">http://www.tobiastiecke.nl/archive/PotassiumProperties.pdf</a> . . . . .	4
<sup>4</sup>	■ model number . . . . .	4
<sup>5</sup>	■ 1:100 or 100:1 . . . . .	4
<sup>6</sup>	■ current? . . . . .	4
<sup>7</sup>	■ how to cite figure . . . . .	5
<sup>8</sup>	■ how to cite figure . . . . .	5
<sup>9</sup>	■ current empty or '-' symbol . . . . .	6
<sup>10</sup>	■ how to cite source . . . . .	6
<sup>11</sup>	■ find name of big yellow probe . . . . .	7
<sup>12</sup>	■ somewhere 2.5-4, find exact value . . . . .	7
<sup>13</sup>	■ Figure: Figure of SHV & BNC connector cable. . . . .	7
<sup>14</sup>	■ 100:1 or 1:100 . . . . .	7
<sup>15</sup>	■ check if really $14\Omega$ or if it event exists . . . . .	10

# References

- [1] Aviv Keshet and Wolfgang Ketterle. "A Distributed, GUI-based, Computer Control System for Atomic Physics Experiments". In: *Review of Scientific Instruments* 84.1 (2013), p. 015105.
- [2] Aviv Keshet. *The Cicero Word Generator*. URL: <https://github.com/akeshet/Cicero-Word-Generator> (visited on 02/20/2020).
- [3] Aviv Keshet. *Cicero Word Generator Technical and User Manual*. URL: <http://akeshet.github.io/Cicero-Word-Generator/Cicero%20Technical%20and%20User%20Manual.pdf> (visited on 02/24/2020).
- [4] Microsoft. *Download .NET (Linux, macOS, and Windows)*. URL: <https://dotnet.microsoft.com/download> (visited on 02/24/2020).
- [5] National Instruments. *NI Driver Downloads - National Instruments*. URL: <https://www.ni.com/en-us/support/downloads/drivers.html> (visited on 02/24/2020).
- [6] National Instruments. *NI Software Gives C++ Runtime Error “Terminated in an Unusual Way” - National Instruments*. URL: <https://knowledge.ni.com/KnowledgeArticleDetails?id=kA00Z0000019YOnSAM&l=en-US> (visited on 02/24/2020).
- [7] Rohde & Schwarz. *HM 507*. URL: [https://cdn.rohde-schwarz.com/hameg-archive/HM507\\_english.pdf](https://cdn.rohde-schwarz.com/hameg-archive/HM507_english.pdf) (visited on 03/28/2020).
- [8] Frank Philipse. *D14363GY123*. URL: <https://frank.pocnet.net/sheets/186/d/D14363GY123.pdf> (visited on 03/10/2020).
- [9] FuG Elektronik GmbH. *HVPS Series HCP*. URL: [https://www.fug-elektronik.de/wp-content/uploads/pdf/Datasheets/EN/HCP\\_data\\_sheet.pdf](https://www.fug-elektronik.de/wp-content/uploads/pdf/Datasheets/EN/HCP_data_sheet.pdf) (visited on 03/23/2020).
- [10] Myrra SAS. *Data Sheet 44231*. URL: <http://www.farnell.com/datasheets/92205.pdf> (visited on 04/07/2020).
- [11] Murata Power Solutions. *ACA-20PC manual*. URL: <https://www.murata.com/products/productdata/8807023804446/aca20pc.pdf> (visited on 04/07/2020).